

Experiment Title: Comparison of Measurement Tools

Objective

1. To compare the accuracy and precision of measurements using a screw gauge and a Vernier caliper by measuring the dimensions of a small object.
2. To compare the accuracy and precision of measurements using a Vernier caliper and a meter scale by measuring the dimensions of a larger object.

Apparatus Required

- Screw Gauge
- Vernier Caliper
- Meter Scale
- Small object (e.g., a small metallic sphere or wire)
- Larger object (e.g., a rod or block)

Theory

1. Screw Gauge:

- Least Count = Pitch/Number of divisions on the circular scale.
- Measures small dimensions like thickness with high precision (usually 0.01 mm).

2. Vernier Caliper:

- Least Count = (Value of one main scale division - Value of one vernier scale division).
- Measures internal, external dimensions, and depth (typically precise to 0.1 mm).

3. Meter Scale:

- Least Count = 1 mm.
- Measures larger dimensions with lower precision.

Procedure

Part A: Comparing Screw Gauge and Vernier Caliper

1. Measure the Length, breadth and height of the small object (e.g., a metallic sphere or block) using the screw gauge:
 - Note the pitch and least count of the screw gauge.

- Record the readings of the main scale and circular scale.
2. Measure the same Length, breadth and height using the Vernier caliper:
 - Note the least count of the Vernier caliper.
 - Record the main scale and Vernier scale readings.

Part B: Comparing Vernier Caliper and Meter Scale

1. Measure the length, breadth or height of the larger object (e.g., a metallic block) using the Vernier caliper:
 - Record the main scale and Vernier scale readings.
2. Measure the same dimensions using the meter scale:
 - Record the readings directly in centimeters or millimeters.

Observations

Screw Gauge Readings

Obs.	MSR	HSR	Corrected HSR	Length (mm) + Corrected HSR × LC)	Average Length (\bar{x}_i)	Deviation ² $(x_i - \bar{x}_i)^2$	$\sigma_{\text{Length}} = \sqrt{(x_i - \bar{x}_i)^2}$
1							
2							
3							
⋮							

Mean Length: $\bar{d} = \dots \pm \sigma_d$

ω

Vernier Caliper Readings

Obs.	MSR	VSR	Length (mm) (MSR + VSR × LC)	Average (\bar{x}_i)	Length	Deviation ² $(x_i - \bar{x}_i)^2$	$\sigma_{\text{Length}} = \sqrt{(x_i - \bar{x}_i)^2}$	=
1								
2								
3								
⋮								

Screw Gauge Readings

Obs.	MSR	HSR	Corrected HSR	Breadth (MSR + HSR × LC)	Average Breadth (\bar{x}_i)	Deviation ² $(x_i - \bar{x}_i)^2$	σ_{Breadth} $\sqrt{\frac{(x_i - \bar{x}_i)^2}{n}}$	=
1								
2								
3								
⋮								

Mean Breadth: $\bar{d} = \dots \pm \sigma_d$

Vernier Caliper Readings

Obs.	MSR	VSR	Breadth (mm) (MSR + VSR × LC)	Average Breadth (\bar{x}_i)	Deviation ² $(x_i - \bar{x}_i)^2$	σ_{Breadth} $\sqrt{\frac{(x_i - \bar{x}_i)^2}{n}}$	=
1							
2							
3							
⋮							

Screw Gauge Readings

Obs.	MSR	HSR	Corrected HSR	Height (mm) + Corrected HSR × LC)	Average Height ($\overline{x_i}$)	Deviation ² $(x_i - \overline{x_i})^2$	$\sigma_{\text{Height}} = \sqrt{(x_i - \overline{x_i})^2}$
1							
2							
3							
⋮							

Mean Height: $\bar{d} = \dots \pm \sigma_d$

Vernier Caliper Readings

Obs.	MSR	VSR	Height (mm) (MSR + VSR × LC)	Average Height ($\overline{x_i}$)	Deviation ² $(x_i - \overline{x_i})^2$	$\sigma_{\text{Height}} = \sqrt{(x_i - \overline{x_i})^2}$
1						
2						
3						
⋮						

0.1 Measurement using meter scale

Obs.	Height (mm)	Average Height (\bar{x}_i)	Deviation ² ($x_i - \bar{x}_i$) ²	$\sigma_{\text{Height}} = \sqrt{(x_i - \bar{x}_i)^2}$
1				
2				
3				
\vdots				

Obs.	Length (mm)	Average Length (\bar{x}_i)	Deviation ² ($x_i - \bar{x}_i$) ²	$\sigma_{\text{Length}} = \sqrt{(x_i - \bar{x}_i)^2}$
1				
2				
3				
\vdots				

Obs.	Breadth (mm)	Average Breadth (\bar{x}_i)	Deviation ² ($x_i - \bar{x}_i$) ²	$\sigma_{\text{Breadth}} = \sqrt{(x_i - \bar{x}_i)^2}$
1				
2				
3				
\vdots				

Calculations

- Least count calculations for the screw gauge and Vernier caliper.
- Error analysis and comparison of precision.

1 Error Analysis

For a set of measurements we estimate the error of the measured value using standard deviation σ which is nothing but

$$\sigma_x = \sqrt{(x_i - \bar{x})^2} \quad (1)$$

1.1 Error associated with measurement using Vernier Calipers

... ..

1.2 Error associated with measurement using Screw Gauge

Result

1. **For small object:** The screw gauge provides higher precision compared to the Vernier caliper for measuring small dimensions.
2. **For larger object:** The Vernier caliper provides better precision compared to the meter scale.

Conclusion

- Screw gauge is more precise for smaller objects, while the Vernier caliper is versatile for medium dimensions.
- Meter scales are less precise and suitable only for approximate measurements of larger objects.

Precautions

- Avoid parallax error while taking readings.
- Ensure zero error correction for both screw gauge and Vernier caliper.
- Handle the instruments carefully to avoid damage.